Scientific correspondence

Free of acrylamide sodium dodecyl sulphate (SDS)-based tissue clearing (FASTClear): a novel protocol of tissue clearing for three-dimensional visualization of human brain tissues

In recent years, advances in laser microscopy and endogenous fluorescent tagging techniques have led to the development of many tissue-clearing strategies, which render tissues optically transparent, allowing large blocks of unsectioned tissue to be visualized in three dimensions (3D). CLARITY (Clear Lipid-Acrylamide-hybridized Rigid Imaging/ exchanged Immunostaining/In situ hybridization-compatible Tissue hYdrogel) is one of the tissue-clearing techniques which works by fixation/hybridization of brain tissue using hydrogel cross-links and subsequent detergentbased delipidation to turn the tissue transparent [1]. Since CLARITY enables molecular probing using immunofluorescence, this technique was deemed suitable for post mortem human brain tissues to demonstrate the potential in visualizing pathologies in Alzheimer's [2], Parkinson's [3] and neurodevelopmental disorders [1,4] in 3D. Early attempts using CLARITY were performed on cortical tissue from the brains of children [1,4], which is considered to be less myelinated than adults, or on thinly sectioned tissues of up to 500 μ m in thickness [1,2]. In our previous study, we have demonstrated that CLARITY can successfully render larger blocks of tissue (about 3 mm in thickness) optically transparent in multiple cortical and subcortical regions of the human brain [3]. However, we noticed the speed of tissue clearing differs between regions depending on the degree of myelination, and duration of formalin fixation. Densely myelinated regions such as the brainstem and spinal cord in archival formalin-fixed tissues could not be rendered transparent with CLARITY. Also we, and other groups, have reported tissue expansion after tissue clearing with CLARITY [5]. Although it was often claimed that the tissue expansion is a transient effect which will be adjusted by subsequent refractive index matching [6], the effects in human brain tissues, especially after

prolonged (>40 days) passive tissue clearing, appeared to be irreversible [3]. In addition, immunolabelling with antibodies, particularly on larger samples, remains challenging because the depth of antibody penetration is still limited [7]. Many current existing tissue-clearing protocols are now available and attempts have been made to improve tissue clearing by combining various protocols, such as CUBIC with RIMS in CLARITY and FRUIT (SeeDB with Scale) [8]. Our aim was to develop an improved and simplified protocol for tissue clearing in the human brain.

The use of acrylamide hydrogel in the CLARITY protocol poses a number of problems including those described above. Tissues embedded with acrylamide hydrogel undergo expansion upon sodium dodecyl sulphate (SDS) clearing and become more fragile as structural integrity is lost [3,5,9]. In addition, when transcardial perfusion cannot be performed, diffusion of hydrogel monomers within a large block of tissue may be limited, leading to incomplete tissue hydrogel hybridization [9]. Furthermore, although pores in polyacrylamide matrices aid lipid exchange and antibody penetration can be enhanced by changing composition of the hydrogel [5], immunolabelling distance appeared to be better in formaldehyde-fixed, unhybridized tissue (without acrylamide hydrogel) compared with those hybridized with both formaldehyde and acrylamide [6]. The original CLARITY study suggested cross-links between hydrogel and formaldehyde aid fixation of protein and nucleic acids during the delipidation process. Protein loss was reported to be significantly lower in acrylamide-embedded tissues compared with unhybridized tissue after SDS clearing [1,5,6]. However, there has been no evidence showing the existence of hydrogel-formaldehyde cross-links. Also, in our recently published study, no significant protein loss was found after SDS clearing of formalin-fixed human brain tissues [10]. Besides, protein loss does not necessarily compromise the quality of immunostaining. Hence, as long as the tissue is well fixed in formaldehyde, we recommend the tissue-clearing procedure to be simplified by omitting the use of acrylamide-based hydrogel.

 $^{\odot}$ 2016 The Authors. Neuropathology and Applied Neurobiology published by John Wiley & Sons Ltd

on behalf of British Neuropathological Society..

This is an open access article under the terms of the Creative Commons Attribution License, which permits use, distribution and reproduction in any medium, provided the original work is properly cited.

Prolonged formaldehyde fixation impedes clearing speed and immunolabelling due to excessive formaldehyde cross-links on the tissue. One simple solution is to use fresh tissue instead of formalin-fixed material. Densely myelinated regions such as brainstem, spinal cord and cortical white matter are difficult to render transparent with SDS delipidation [3]. Reiner and colleagues introduced iDISCO which demonstrated that it is possible to immunostain a piece of optically opaque, formaldehyde-fixed tissue by permeabilizing it in a cocktail of detergent [11]. The tissue can then be rendered optically transparent with organic solvents using the 3DISCO clearing technique [12]. Apart from its delipidation properties, SDS is known to be a protein denaturant, which accounts for its use for antigen retrieval in traditional immunohistochemistry. As a result, to improve immunolabelling in tissue clearing, the combination of SDS delipidation in CLARITY and detergent permeabilization in iDISCO can be useful. In addition, we observed that for densely expressed antigens, such as glial fibrillary acidic protein GFAP and neurofilament, the use of a low antibody concentration (1:1000) with daily supplement (to a concentration of 1:100-1:50) may be useful to prevent antibodies being 'trapped' at tissue surface.

With the omission of acrylamide-hydrogel in CLARITY and combination with the iDISCO technique, we introduce FASTClear (Free of Acrylamide SDS-based Tissue Clearing) for the immunostaining and threedimensional visualization of human brain tissue (Table 1). First, fresh human brain blocks of up to 1 cm in thickness were fixed in 4% paraformaldehyde or 10% neutral-buffered formalin at 4°C for 3 days. Then, the tissue was trimmed to about 3 mm in thickness (the maximum working distance of confocal objectives) before delipidation in 4% SDS-boric acid buffer at 50°C for a minimum of 5 days. Best immunostaining was achieved if the region of interest is rendered optically transparent at this stage. However, it is possible to proceed to immunostaining even if the tissues do not achieve full transparency at this point. The tissue was then washed thoroughly in phosphate-buffered saline (PBS) with 0.1% Triton-X 100 (PBS-Triton) at 50°C $(3 \times 1 h)$. Then, the tissue was permeabilized and blocked in 0.6 M glycine, 0.2% Triton X-100, 6% Donkey serum, 20% dimethyl sulfoxide (DMSO) dissolved in PBS overnight at 37°C. Next, after washing the tissue in PBS-Triton for 2×1 h at 37°C, it was incubated with primary antibody diluted in 0.2% Tween-20, 5% DMSO, 3% Donkey serum, 0.01% sodium azide in PBS for a minimum of 2 days at 37°C. Following another wash in PBS-Triton $(3 \times 1 h, then overnight incuba$ tion at 37°C), the tissue was incubated with a secondary antibody conjugated with Alexa Fluor[©] fluorophores diluted in the same diluent as above for the same number of days as primary antibody incubation at 37°C. A nuclear counterstain, 4',6-diamidino-2phenylindole (1:100 from a stock of 1 µg/ml diluted with 1:1 water: DMSO), can be added at this stage. After that, the tissue was washed thoroughly in PBS-Triton $(5 \times 1 \text{ hr}; \text{ then overnight})$ and proceeded to refractive index matching. For tissues that have been rendered optically transparent at the delipidation step. immersion in 47% 2,2'-thiodiethanol (vol/vol) in 0.01 M PBS without saline or 70% w/v sorbitol in 0.1 M phosphate buffer as previously described [3] was done for refractive index matching. For tissues that did not achieve transparency at the delipidation step (or if microscope objectives are designed for high refractive index solution), they can be dehydrated and refractive index matched as per the 3DISCO protocol [12]. Briefly, tissue was immersed in 50% tetrahydrofuran (THF) (overnight), 70% THF (1 h), 80% THF (1 h), 100% THF (1 h), 100% THF (1 h) and finally dibenzyl ether until the tissue became optically transparent. Tissues were then mounted and visualized using a singlephoton or two-photon confocal microscope.

We applied this protocol on a piece of fresh spinal cord tissue, which is difficult to render transparent with the traditional CLARITY technique, and successfully visualized the three-dimensional structure of the ventral root to a depth of 508.52 µm using immunostaining for neurofilament (Figure 1a). However, it has to be noted that there is heterogeneity in the penetration depth of the antibody with some parts of the spinal cord with labelling to 89 µm only. The reasons behind the variation in penetration depth within the same piece of tissue remain largely unknown, and we hypothesize that local tissue structural biochemical properties or vasculature can affect antibody penetration. Furthermore, we have also demonstrated the potential of FASTClear with other antibodies such as tyrosine hydroxylase (Millipore AB152, UK) and microtubule-associated protein 2 (abcam ab5392, UK) (data not shown) and on formalinfixed tissue in a tissue bank for the 3D visualization of Purkinje neurons in the cerebellum (Figure 1b), with a depth of immunolabelling to 66.5 µm.

© 2016 The Authors. Neuropathology and Applied Neurobiology published by John Wiley & Sons Ltd on behalf of British Neuropathological Society.

NAN 2017; 43: 346-351

Table	1. Free of	acrylamide	sodium	dodecyl	sulphate	(SDS)-based	tissue-clearing	(FASTClear)	protocol
-------	------------	------------	--------	---------	----------	-------------	-----------------	-------------	----------

Step	Time
 Fixation in 10% neutral-buffered formalin/4% paraformaldehyde (PFA) N.B. fixation time depends on size of fresh tissue block. Typically, a 1-cm-thick block will take around 3 days to be fully fixed. This step is only required for fresh tissue. Proceed to Step 2 for formalin-fixed tissue. Caution: use a tightly sealed container (±parafilm) for fixation as formalin/PFA is toxic. 	Minimum of 3 days at 4°C
2. Dissect into smaller blockA maximum of 3 mm in thickness is recommended due to immunolabelling diffusion and confocal objectives working distance limits. Note that the sectioning surface should be as flat as possible and designed to be the future imaging surface.	
 3. Immerse in 4% SDS buffer This step improves antibody labelling and it is recommended the tissue is immersed in SDS buffer until transparency is reached. For prolonged fixed tissue (>2 years in fixation), a 2-mm-thick block can reach transparency in 3 months. [However, do note that complete transparency of tissue is not a necessity as tissue will become transparent at the final refractive index matching step] Frequent change in buffer (daily to twice weekly) can improve the speed of reaching transparency. 	Minimum of 5 days at 50°C oven
4. Washing in 0.1% phosphate-buffered saline (PBS)-Triton	3×1 h at 50°C
 5. Blocking and permeabilization in blocking medium [0.6 M glycine, 0.2% Triton X-100, 6% Donkey serum, 20% dimethyl sulfoxide (DMSO) in PBS] Add enough blocking medium to cover the tissue. Optional if the antibody is known to be of high specificity. 	Overnight at 37°C
6. Washing in 0.1% PBS-Triton	2 x 1 h at 37°C
 7. Primary antibody incubation (diluted in 0.2% Tween-20, 5% DMSO, 3% Donkey serum, 0.01% sodium azide in PBS) Start with a low concentration (e.g. 1:1000; 2 μl in 2 ml of diluent), supplementing antibody daily/twice daily until a final concentration of around 1:50–1:100 is reached. Optimal concentration and days of incubation vary between antibodies. As an example, tyrosine hydroxylase antibodies (Millipore AB152) can reach complete penetration to a depth of 2 mm on each side in 3 days at a final concentration of 1:100. If multiple antigen labelling is required, it is recommended to perform immunolabelling sequentially. 	Minimum of 2 days at 37°C
8. Washing in 0.1% PBS-Triton	3×1 h at 37°C; then overnight at 37°C
 9. Secondary antibody incubation (diluted in 0.2% Tween-20, 5% DMSO, 3% Donkey serum, 0.01% sodium azide in PBS) • Same as Step 7 above. • 4',6-Diamidino-2-phenylindole or fluorophore-conjugated lectin can be added at this stage (1:100 from a stock of 1 μg/ml diluted with 1:1 water: DMSO) for better tissue orientation. 	Minimum of 2 days at 37°C
10. Washing in 0.1% PBS-Triton	5×1 h at 37°C; then overnight at 37°C
 11. Immersion in refractive index matching medium If tissue is transparent/almost transparent after Step 3 Immerse in 47% 2,2'-thiodiethanol diluted in 0.01 M PBS without saline or 70% w/v Sorbitol in 0.1 M phosphate buffer (as previously described in ref. 3). If tissue is opaque after Step 3, follow 3DISCO clearing method: Dehydrate tissue in 50% tetrahydrofuran (THF) (overnight), 70% THF (1 h), 80% THF (1 h), 100% THF (1 h) 100% THF (1 h) then dibenzyl ether until transparency is reached. Caution: perform this step in the fume bood. 	RT until transparency is reached

FASTClear is a greatly simplified and more userfriendly tissue-clearing protocol for human brain tissues which reduces the overall processing time from tissue fixation to immunostaining and visualization to a minimum of 16 days (Figure 2). However, it is still not possible to achieve immunolabelling to the full thickness

© 2016 The Authors. Neuropathology and Applied Neurobiology published by John Wiley & Sons Ltd NAN 2 on behalf of British Neuropathological Society.



Figure 1. Human brain tissue processed and immunostained using free of acrylamide sodium dodecyl sulphate (SDS)-based tissue clearing (FASTClear). (a) z-Stack image of a ventral root of a piece of spinal cord fresh tissue (3 mm thick) immunostained using antineurofilament primary antibody (final concentration 1:100; Dako M0762) and Alexa-fluor 488 conjugated donkey anti-mouse secondary antibody. Stained tissue was visualized using a Zeiss 780 inverted confocal microscope (Carl Zeiss, Germany) with ×10 objective (imaging depth to 508.519 μ m, z-stack step size 3.03 μ m). (b) z-Stack image of a piece of fixed cerebellar tissue immunostained using antibodies against neurofilament (green; final concentration 1:100; Dako M0762) and βIII-tubulin (red; final concentration 1:100; Millipore AB9354) and counterstained with 4',6-diamidino-2-phenylindole (blue). Stained tissue was visualized using a Leica SP5 (Leica Microsystems, UK) confocal microscope with ×40 objective (imaging depth to 66.5 μ m, z-stack step size 0.38 μ m).



Figure 2. Workflow diagram of free of acrylamide sodium dodecyl sulphate (SDS)-based tissue clearing (FASTClear).

of the tissue and this technique is yet to be optimized for archival formalin-fixed tissues. Also, tissue that requires clearing with organic solvent may undergo shrinkage. Although fine structures are likely to be preserved, a high-power objective will be required which often has a lower working distance. Nevertheless, with the development of smaller molecular probes such as nanobodies and aptamers [7], further improvement of human brain tissue clearing will be made, leading to a new era of three-dimensional histology.

© 2016 The Authors. Neuropathology and Applied Neurobiology published by John Wiley & Sons Ltd on behalf of British Neuropathological Society.

Acknowledgements

The work was funded through an innovation grant (K-1403) from Parkinson's UK and the Alzheimer's Research UK Imperial College London Network Centre Grant (ARUK-2014NC-IMP). The authors would like to thank Mr Steve Rothery and other staff at the Imperial FILM facility, and technical staff of the Parkinson's UK Tissue Bank for their assistance. Finally, we express our deepest appreciation to the donors and their families for donating human brain tissue for research.

Author contributions

All authors contributed to the experimental design and conceived the study. A.K.L.L. and H.M.L. carried out all the experimental work, drafted and revised the manuscript. R.C.C.C. and S.M.G. supervised the research. All authors read, reviewed and edited the final manuscript.

Ethical considerations

The work conducted on human tissue was under ethical approval held by the Parkinson's UK Brain Bank at Imperial College London [Registered charity in England and Wales (258197) and in Scotland (SC037554); Multicentre Research Ethics Committee approval reference number: 07/MRE09/72]. Parkinson's UK Brain Bank is an approved Research Tissue Bank by the Wales Research Ethics Committee (Ref. No. 08/MRE09/ 31+5). Informed consent was obtained prospectively for the use of *post mortem* brain tissues and brain samples were obtained and prepared in accordance to the Wales Research Ethics Committee approved protocols.

Conflict of interest

All authors declare no conflict of interest.

A. K. L. Liu*,⁺¹ H. M. Lai*,^{±1} R. C.-C. Chang⁺,[±], S. M. Gentleman* *Neuropathology Unit, Division of Brain Sciences, Department of Medicine, Imperial College London, London, UK, ⁺Laboratory of Neurodegenerative Diseases, LKS Faculty of Medicine,

School of Biomedical Sciences, ‡State Key Laboratory of Brain and Cognitive Sciences, and §LKS Faculty of Medicine, Research Centre of Heart, Brain, Hormone, and Healthy Aging, The University of Hong Kong, Pokfulam, Hong Kong SAR E-mail: king.liu09@imperial.ac.uk ¹These authors contributed equally to this work.

References

- 1 Chung K, Wallace J, Kim S-Y, Kalyanasundaram S, Andalman AS, Davidson TJ, Mirzabekov JJ, Zalocusky KA, Mattis J, Denisin AK, Pak S, Bernstein H, Ramakrishnan C, Grosenick L, Gradinaru V, Deisseroth K. Structural and molecular interrogation of intact biological systems. *Nature* 2013; **497**: 332–7
- 2 Ando K, Laborde Q, Lazar A, Godefroy D, Youssef I, Amar M, Pooler A, Potier M-C, Delatour B, Duyckaerts C. Inside Alzheimer brain with CLARITY: senile plaques, neurofibrillary tangles and axons in 3-D. *Acta Neuropathol* 2014; **29**: 457–9
- 3 Liu AKL, Hurry ME, Ng OT-W, DeFelice J, Lai HM, Pearce RK, Wong GT-C, Chang RC-C, Gentleman SM. Bringing CLARITY to the human brain: visualisation of Lewy pathology in three-dimensions. *Neuropathol Appl Neurobiol* 2015; **42**(6): 573–87.
- 4 Costantini I, Ghobril J-P, Di Giovanna AP, Mascaro ALA, Silvestri L, Müllenbroich MC, Onofri L, Conti V, Vanzi F, Sacconi L, Guerrini R, Markram H, Iannello G, Pavone FS. A versatile clearing agent for multimodal brain imaging. *Sci Rep* 2015; **5**: 9808
- 5 Yang B, Treweek JB, Kulkarni RP, Deverman BE, Chen CK, Lubeck E, Shah S, Cai L, Gradinaru V. Single-cell phenotyping within transparent intact tissue through whole-body clearing. *Cell* 2014; **158**: 945– 58
- 6 Treweek JB, Chan KY, Flytzanis NC, Yang B, Deverman BE, Greenbaum A, Lignell A, Xiao C, Cai L, Ladinsky MS, Bjorkman PJ, Fowlkes CC, Gradinaru V. Whole-body tissue stabilization and selective extractions via tissue-hydrogel hybrids for high-resolution intact circuit mapping and phenotyping. *Nat Protoc* 2015; **10**: 1860–96
- 7 Marx V. Optimizing probes to image cleared tissue. *Nat Methods* 2016; **13**: 205–9
- 8 Susaki EA, Ueda HR. Whole-body and whole-organ clearing and imaging techniques with single-cell resolution: toward organism-level systems biology in mammals. *Cell Chem Biol* 2016; **23**: 137–57
- 9 Murray E, Cho JH, Goodwin D, Ku T, Swaney J, Kim SY, Choi H, Park YG, Park JY, Hubbert A, McCue M, Vassallo S, Bakh N, Frosch MP, Wedeen VJ, Seung HS, Chung K. Simple, scalable proteomic imaging for highdimensional profiling of intact systems. *Cell* 2015; 163: 1500–14
- 10 Lai HM, Liu AKL, Ng W-L, DeFelice J, Lee WS, Li H, Li W, Ng HM, Chang RC-C, Lin B, Wu W, Gentleman

© 2016 The Authors. Neuropathology and Applied Neurobiology published by John Wiley & Sons Ltd on behalf of British Neuropathological Society.

SM. Rationalisation and validation of an acrylamide-free procedure in three-dimensional histological imaging. *PLoS ONE* 2016; 11: e0158628.

- 11 Renier N, Wu Z, Simon DJ, Yang J, Ariel P, Tessier-Lavigne M. iDISCO: a simple, rapid method to immunolabel large tissue samples for volume imaging. *Cell* 2014;**159**:896–910.
- 12 Ertürk A, Becker K, Jährling N, Mauch CP, Hojer CD, Egen JG, Hellal F, Bradke F, Sheng M, Dodt H-U.

Three-dimensional imaging of solvent-cleared organs using 3DISCO. *Nat Protoc* 2012; 7: 1983–95

Received 20 May 2016 Accepted after revision 25 August 2016 Published online Article Accepted on 7 September 2016