



# Whole-Genome Sequence of Staphylococcus epidermidis Tü3298

### Josephine C. Moran, Malcolm J. Horsburgh

Institute of Integrative Biology, University of Liverpool, Liverpool, United Kingdom

*Staphylococcus epidermidis* Tü3298 is a frequently used laboratory strain, known for its production of epidermin and absence of the *icaABCD* operon. We report the whole-genome sequence of this strain, a 2.5-kb genome containing 2,332 genes.

Received 26 January 2016 Accepted 29 January 2016 Published 10 March 2016

**Citation** Moran JC, Horsburgh MJ. 2016. Whole-genome sequence of *Staphylococcus epidermidis* Tü3298. Genome Announc 4(2):e00112-16. doi:10.1128/genomeA.00112-16.

**Copyright** © 2016 Moran and Horsburgh. This is an open-access article distributed under the terms of the Creative Commons Attribution 4.0 International license. Address correspondence to Malcolm J. Horsburgh, M.J.Horsburgh@liverpool.ac.uk.

**S***taphylococcus epidermidis* is a frequent cause of hospitalacquired infections, exhibiting increasing levels of antimicrobial resistance that are a growing concern. The bacterium primarily causes infections of catheters and other inserted medical devices. The life cycle of S. epidermidis as a normal and dominant resident of the human skin flora affords this opportunistic pathogen with access to these inserted medical devices that is difficult to prevent.

As a colonizing species of human skin, S. epidermidis must compete with other bacteria for space and nutrients (1, 2). Many staphylococci produce lantibiotics, antibiotic-like peptides, as a competition mechanism (3). S. epidermidis Tü3298 produces the lantibiotic epidermin, active against many Grampositive bacteria (4). This was the first lantibiotic shown to be post-translationally modified, and remains one of the bestcharacterized staphylococcal lantibiotics, with its biosynthesis, regulation, mode of action, and the nature of bacterial resistance all having been elucidated (3). The optimized protocol for transformation by electroporation of S. epidermidis Tü3298 (5) made this strain ideal for use in studies of gene function via mutagenesis. Subsequently, S. epidermidis Tü3298 has been used in studies of the role of agr (6–8),  $\sigma^{B}$  (9), and SepA (10). This strain is also widely used as an *icaABCD* negative control in biofilm assays (11, 12).

Libraries were prepared for sequencing with Nextera DNA kits (Illumina) and were sequenced on the Illumina MiSeq platform. Contigs were assembled using VelvetOptimiser (Velvet version 1.2.06) and were annotated using Prokka (13). Sequences were assembled into 162 contigs; the average coverage for the assembly was 14.

The *S. epidermidis* Tü3298 genome is 2,459,658 bp. It contains 2,332 protein coding sequences, 360 of these are hypothetical.

The genome sequence of *S. epidermidis* Tü3298 will contribute to easier genetic manipulation of this strain, and will enable further sequence-dependent studies in the future.

**Nucleotide sequence accession numbers.** This draft genome sequence has been deposited in the ENA under the accession no. CZRO00000000. The version described in this paper is the first version, CZRO02000000.

## ACKNOWLEDGMENT

DNA was sequenced at the Center for Genomic Research, Liverpool.

## FUNDING INFORMATION

This work, including the efforts of Josephine Moran, was funded by Biotechnology and Biological Sciences Research Council (BBSRC) (BB/ I532161/1 and BB/L023040/1).

The funders were not involved in the study design, collection of samples, analysis of data, interpretation of data or the writing of this report or the decision to submit this report for publication.

## REFERENCES

- Coates R, Moran J, Horsburgh MJ. 2014. Staphylococci: colonizers and pathogens of human skin. Future Microbiol 9:75–91. http://dx.doi.org/ 10.2217/fmb.13.145.
- Cogen AL, Nizet V, Gallo RL. 2008. Skin microbiota: a source of disease or defence? Br J Dermatol 158:442–455. http://dx.doi.org/10.1111/j.1365 -2133.2008.08437.x.
- 3. Götz F, Perconti S, Popella P, Werner R, Schlag M. 2014. Epidermin and gallidermin: staphylococcal lantibiotics. Int J Med Microbiol 304:63–71. http://dx.doi.org/10.1016/j.ijmm.2013.08.012.
- Augustin J, Rosenstein R, Wieland B, Schneider U, Schnell N, Engelke G, Entian KD, Götz F. 1992. Genetic analysis of epidermin biosynthetic genes and epidermin-negative mutants of *Staphylococcus epidermidis*. Eur J Biochem 204:1149–1154. http://dx.doi.org/10.1111/j.1432 -1033.1992.tb16740.x.
- Augustin J, Götz F. 1990. Transformation of *Staphylococcus epidermidis* and other staphylococcal species with plasmid DNA by electroporation. FEMS Microbiol Lett 54:203–207.
- Vuong C, Gerke C, Somerville GA, Fischer ER, Otto M. 2003. Quorumsensing control of biofilm factors in *Staphylococcus epidermidis*. J Infect Dis 188:706–718. http://dx.doi.org/10.1086/377239.
- Vuong C, Götz F, Otto M. 2000. Construction and characterization of an agr deletion mutant of *Staphylococcus epidermidis*. Infect Immun 68: 1048–1053. http://dx.doi.org/10.1128/IAI.68.3.1048-1053.2000.
- Otto M, Süssmuth R, Jung G, Götz F. 1998. Structure of the pheromone peptide of the *Staphylococcus epidermidis* agr system. FEBS Lett 424: 89–94. http://dx.doi.org/10.1016/S0014-5793(98)00145-8.
- Kies S, Otto M, Vuong C, Götz F. 2001. Identification of the sigB operon in *Staphylococcus epidermidis*: construction and characterization of a sigB deletion mutant. Infect Immun 69:7933–7936. http://dx.doi.org/10.1128/ IAI.69.12.7933-7936.2001.
- Teufel P, Götz F. 1993. Characterization of an extracellular metalloprotease with elastase activity from *Staphylococcus epidermidis*. J Bacteriol 175:4218–4224.

- Tormo MA, Knecht E, Götz F, Lasa I, Penadés JR. 2005. Bap-dependent biofilm formation by pathogenic species of *Staphylococcus*: evidence of horizontal gene transfer? Microbiology 151:2465–2475. http://dx.doi.org/ 10.1099/mic.0.27865-0.
- Handke LD, Slater SR, Conlon KM, O'Donnell ST, Olson ME, Bryant KA, Rupp ME, O'Gara JP, Fey PD. 2007. SigmaB and SarA indepen-

dently regulate polysaccharide intercellular adhesin production in *Staphylococcus epidermidis*. Can J Microbiol 53:82–91. http://dx.doi.org/ 10.1139/w06-108.

 Seemann T. 2014. Prokka: rapid prokaryotic genome annotation. BioInformatics 30:2068-2069. http://dx.doi.org/10.1093/ bioinformatics/btu153.